

### Clinical Equipment Risk Evaluation Form

Equipment Type Description:

Manufacturer:

Environment Area Used:

#### Evaluation Information:

Risk Category I: Equipment Function (E) Score:

Risk Category II: Clinical Application (A) Score:

Risk Category III: PM Requirement (P) Score:

Risk Category IV: Environmental Use Area: (U) Score:

#### Evaluation Scoring:

Total =  $E + A + [(P + U) / 2]$  Score:

#### Inventory Classification Results:

- Priority I (Score is in the 18 - 20 range)
- Priority II (Score is in the 15 - 17 range)
- Priority III (Score is in the 12 - 14 range)
- Hazard surveillance check only (score is in the 6 - 11 range)
- Deleted from the management program inventory (Score is 5 or less)

#### Additional Comments:

#### Approval and Responsibility

Reviewed by:  Date:

Approved by:  Date: